



Anti-Cytokeratin polyclonal antibody (DPATB-H83439)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-wide spectrum Cytokeratin Polyclonal Antibody
Target	wide spectrum Cytokeratin
Immunogen	Cytokeratin isolated from bovine muzzle epidermis.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr, ICC/IF, ICC
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.5.
Preservative	None
Storage	Store at +4°C.

BACKGROUND

Introduction

Cytokeratins are proteins of keratin-containing intermediate filaments found in the intracytoplasmic cytoskeleton of epithelial tissue. The term "cytokeratin" began to be used in the late 1970s (for example, see "Intermediate-sized filaments of human endo

Keywords

pan cytokeratin,
