



Anti-DsbA L (full length) polyclonal antibody (DPATB-H83440)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Product Overview | Rabbit Anti-DsbA L Polyclonal Antibody |
| Target | DsbA L |
| Immunogen | Recombinant full length Mouse DsbA L protein. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Mouse |
| Purification | Whole antiserum |
| Conjugate | Unconjugated |
| Applications | WB, IHC-Fr |
| Positive Control | 3T3-L1 adipocytes and Mouse white and brown adipose tissue lysate. |
| Format | Liquid |
| Size | 100 µl |
| Buffer | 40% Glycerol, Whole serum |
| Preservative | None |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Disulfide-bond-A oxidoreductase-like protein (DsbA-L, previously named as GST Kappa) is an adiponectin-interacting protein. DsbA-L is highly expressed in adipose tissue, and its expression level is negatively correlated with obesity in mice and humans. Ds |
| Keywords | DSBA like thioredoxin domain containing protein; |