



Anti-PGCF (C-terminal) polyclonal antibody (DPATB-H83448)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PGCF Polyclonal Antibody
Target	PGCF
Immunogen	Synthetic peptide for particulate guanylyl cyclase F. All peptides were amidated for conjugation. The PGCF selective antibodies were prepared against a unique peptide selective for only PGCF that maps at or close to the C terminal region of the protein. (
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Glycerol, Sodium chloride, HEPES
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Cyclic GMP (cGMP), a key messenger in several signal transduction pathways, the intracellular levels of cGMP ate maintained by the activity of opposing enzymes: synthesizing gualylyl cyclases (GC) and hydrolyzing phosphodiesterases (PDEs). The synthesizin
Keywords	Particulate Guanylate Cyclase F;