



Anti-Endomorphin 1 polyclonal antibody (DPATB-H83452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Endomorphin 1 Polyclonal Antibody
Target	Endomorphin 1
Immunogen	Synthetic peptide conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-Fr, ELISA
Format	Liquid
Size	50 µg
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The endogenous agonists for the mu-opiate receptors, endomorphin-1; Tyr-Pro-Trp-Phe-NH₂,

and endomorphin-2; Tyr-Pro-Phe-Phe-NH₂ have been reported to show hypotensive activity and Ca²⁺ channel current inhibitory activity [Biochem. Biophys. Res. Commun.,23

Keywords

EM 1;
