



# Anti-Cytomegalovirus type AD169 Polyclonal antibody (DPATB-H83461)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Goat Anti-Cytomegalovirus AD169 Polyclonal Antibody
Target	Cytomegalovirus type AD169
Immunogen	Purified virions of Cytomegalovirus, Strain AD 169.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Cytomegalovirus
Purification	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, Conjugation, ICC, WB, ICC, IF
Format	Liquid
Size	1 mg
Buffer	PBS.
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Cytomegalovirus (CMV) is a member of the Herpes virus family. Members of this family have a characteristic virion structure. The double stranded DNA genome is contained within an icosahedral capsid which is embedded in a proteinaceous layer (tegument) and

---

**Keywords**

HCMV; HHV5; Human Herpesvirus 5; UL55 (strain AD169);

---