



# Anti-RSV Polyclonal antibody (DPATB-H83469)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                    |
|------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>      | Goat Anti-Respiratory Syncytial Virus Polyclonal Antibody                                                                          |
| <b>Target</b>                | RSV                                                                                                                                |
| <b>Immunogen</b>             | Human RSV isolate                                                                                                                  |
| <b>Isotype</b>               | IgG                                                                                                                                |
| <b>Source/Host</b>           | Goat                                                                                                                               |
| <b>Species Reactivity</b>    | RSV                                                                                                                                |
| <b>Purification</b>          | IgG fraction                                                                                                                       |
| <b>Conjugate</b>             | Unconjugated                                                                                                                       |
| <b>Applications</b>          | ICC/IF, IHC-P                                                                                                                      |
| <b>Cellular Localization</b> | Virion. Host cytoplasm                                                                                                             |
| <b>Format</b>                | Liquid                                                                                                                             |
| <b>Size</b>                  | 100 µl                                                                                                                             |
| <b>Buffer</b>                | 10mg/ml BSA, 0.01M PBS, pH 7.2                                                                                                     |
| <b>Preservative</b>          | 0.1% Sodium Azide                                                                                                                  |
| <b>Storage</b>               | Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Store in the dark. Avoid repeated freeze / thaw cycles. |

## BACKGROUND

**Introduction**

Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en

---

**Keywords**

RS virus,RSV,

---