



Anti-Vaccinia virus Polyclonal antibody (DPATB-H83470)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vaccinia Virus Polyclonal Antibody
Target	Vaccinia virus
Immunogen	Native virus, Lister strain.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Vaccinia virus
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr, ICC/IF
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Vaccinia virus is an Orthopoxvirus, containing double stranded DNA. Fusion protein plays an important role in the entry of enveloped virus into cells. As vaccinia virus has a wide host range,
--------------	---

it is conceivable that certain cellular components that are ubi

Keywords

A27L,Orthopoxvirus,
