



Anti-EBV 3B Polyclonal antibody (DPATB-H83473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-EBNA 3B Polyclonal Antibody
Target	EBV 3B
Immunogen	Full length protein.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	EBV
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	250 µg
Buffer	PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	EBNA 3B is one of only 20 viral proteins expressed in Epstein Barr Virus (EBV) transformed B lymphoblasts cultured in vitro, though the expression of EBNA 3B is not required for EBV-
--------------	--

mediated transformation in vitro. This 165kDa nuclear protein functions as

Keywords

EBNA 4; EBV; EBV nuclear antigen 4; Epstein Barr nuclear antigen 4; HHV4;
