



Magic™ Anti-RSV Polyclonal antibody (DPATB-H83488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Respiratory Syncytial Virus Polyclonal Antibody
Target	RSV
Immunogen	Leukosis RSV-B.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	RSV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Neut
Cellular Localization	Virion. Host cytoplasm
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en
Keywords	RS virus,RSV,