



Anti-IBV Polyclonal antibody (DPATB-H83489)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Guinea pig Anti-Influenza B Virus Polyclonal Antibody
Target	IBV
Immunogen	Intact, heat killed Haemophilus influenza type B.
Isotype	IgG
Source/Host	Guinea pig
Species Reactivity	IBV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Nuclear
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Influenza A and B are the two types of influenza viruses that cause epidemic human disease. Influenza B viruses currently circulating are divided into two antigenically and genetically distinct groups: the Victoria and Yamagata lineages.

Keywords

Influenza virus type B; NP; Nucleocapsid protein; Nucleoprotein; Protein N;
