



# Anti-HSV type 2 Polyclonal antibody (DPATB-H83493)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Guinea pig Anti-HSV2 Polyclonal Antibody
Target	HSV type 2
Immunogen	Cell preparation: HSV type3whole virus
Isotype	IgG
Source/Host	Guinea pig
Species Reactivity	HSV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, Neut
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Herpes simplex type3(HSV2) belongs to a family that includes HSV1, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus. HSV1 and HS10 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral p??

---

**Keywords**

Herpes Simplex Virus Type 2;

---