



# Anti-Cytomegalovirus Glycoprotein B Polyclonal antibody (DPATB-H83496)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Goat Anti-Cytomegalovirus Glycoprotein B Polyclonal Antibody
Target	Cytomegalovirus Glycoprotein B
Immunogen	cell preparation: Purified Cytomegalovirus virions of Strain AD169
Isotype	IgG
Source/Host	Goat
Species Reactivity	Cytomegalovirus
Conjugate	Unconjugated
Applications	ICC, IF, WB
Cellular Localization	Virus envelope.
Format	Liquid
Size	1 ml
Buffer	10mM PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Human cytomegalovirus (HCMV) glycoprotein B (gB) is an abundant virion envelope protein that has been shown to be essential for the infectivity of HCMV. HCMV gB is also one of the most immunogenic virus-encoded proteins, and a significant fraction of viru

**Keywords**

CMV; CMV glycoprotein B; CMV glycoprotein GP55; Cytomegalovirus glycoprotein GP55; gB; Glycoprotein B; Glycoprotein GP55; HCMV; HCMV glycoprotein B; HCMV glycoprotein GP55; HHV 5; HHV6glycoprotein B; HHV6glycoprotein GP55; HHV5 glycoprotein B; HHV5 glycop