



Anti-Cytomegalovirus Glycoprotein B Polyclonal antibody (DPATB-H83501)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Cytomegalovirus Glycoprotein B Polyclonal Antibody
Target	Cytomegalovirus Glycoprotein B
Immunogen	Cell preparation: purified virions of strain AD169.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Cytomegalovirus
Conjugate	Unconjugated
Applications	ICC, IF, ELISA, WB, IHC-Fr
Cellular Localization	Virus envelope.
Format	Liquid
Size	1 ml
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Human cytomegalovirus (HCMV) glycoprotein B (gB) is an abundant virion envelope protein that has been shown to be essential for the infectivity of HCMV. HCMV gB is also one of the most immunogenic virus-encoded proteins, and a significant fraction of viru

Keywords

CMV; CMV glycoprotein B; CMV glycoprotein GP55; Cytomegalovirus glycoprotein GP55; gB; Glycoprotein B; Glycoprotein GP55; HCMV; HCMV glycoprotein B; HCMV glycoprotein GP55; HHV6glycoprotein B; HHV6glycoprotein GP55; HHV5 glycoprotein B; HHV5 glycoprotein
