



Anti-HCV type 5 NS5 Polyclonal antibody (DPATB-H83509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Hepatitis C Virus genotype6NS5 Polyclonal Antibody
Target	HCV type 5 NS5
Immunogen	HCV NS5a
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HCV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Endoplasmic reticulum membrane
Format	Liquid
Size	100 µg
Buffer	1X PBS
Preservative	0.05% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Hepatitis C Virus is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of about 3,000 amino acids. The polyprotein is processed by host cell and viral protease

Keywords

HCV NS5; HCV NS5A; Hepatitis C virus NS5A; Non structural protein 5A; p56;
