



Anti- US3 Polyclonal antibody (DPATB-H83516)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-US3 Polyclonal Antibody
Target	US3
Immunogen	Synthetic peptide derived from an internal domain of Human herpes simplex virus 1 (HHV1 strain 17) US3 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV 1
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Tegument protein; phosphorylates nuclear egress lamina protein; mediates phosphorylation of HDAC1 and HDAC2 and other cellular and viral proteins.
Keywords	Serine/threonine protein kinase US3,