



Anti- RL2 Polyclonal antibody (DPATB-H83517)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-RL2 Polyclonal Antibody
Target	RL2
Immunogen	A synthetic peptide corresponding to C terminal residues of Human herpesvirus3RL2 protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

BACKGROUND

Introduction	RL2 (Trans acting transcriptional protein ICP0) contains 1 RING type zinc finger. RL2 disrupts ND10 and proteasome dependent degradation of several cellular proteins.
---------------------	--

Keywords

Trans acting transcriptional protein ICP0,VMW118 protein,
