



Anti-RSV Major Glycoprotein Polyclonal antibody (DPATB-H83521)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Respiratory Syncytial Virus Major Glycoprotein Polyclonal Antibody
Target	RSV Major Glycoprotein
Immunogen	Synthetic peptide: LCKSICKTIPSNKPKKKP, corresponding to amino acids 187-198 of Respiratory Syncytial Virus Major Glycoprotein B1 strain. LCKSICKTIPSNKPKKKP Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	RSV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell surface of infected cells.
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en
Keywords	Attachment glycoprotein G,Fusion glycoprotein F1,Major surface glycoprotein G,Membrane-bound glycoprotein,mG,RSV,