



Anti-RSV Polyclonal antibody (DPATB-H83522)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Respiratory Syncytial Virus Polyclonal Antibody
Target	RSV
Immunogen	Respiratory Syncytial Virus isolate, confirmed.
Isotype	IgG
Source/Host	Goat
Species Reactivity	RSV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, ICC/IF
Cellular Localization	Virion. Host cytoplasm
Format	Liquid
Size	500 µg
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en

Keywords

RS virus,RSV,
