



Anti-IAV Polyclonal antibody (DPATB-H83526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Influenza A Virus Hemagglutinin H1 Polyclonal Antibody
Target	IAV
Immunogen	A synthetic peptide from Influenza A Virus Hemagglutinin H1 (Genbank ACA28844).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	IAV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB
Positive Control	Influenza A Virus Hemagglutinin H1
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Influenza A virus is a major public health threat. Novel influenza virus strains caused by genetic drift and viral recombination emerge periodically to which humans have little or no immunity, resulting in devastating pandemics. Influenza A can exist in a

Keywords

HA; Hemagglutinin; Hemagglutinin HA;
