



Magic™ Anti-IAV Polyclonal antibody (DPATB-H83527)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Influenza A Virus Polyclonal Antibody
Target	IAV
Immunogen	Tissue / cell preparation: Influenza A/Puerto Rico 8/34 (H1N1) virus.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	IAV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Nuclear
Format	Liquid
Size	200 µl
Buffer	40% Glycerol, Whole serum.
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Influenza A and B are the two types of influenza viruses that cause epidemic human disease. Influenza type C infections cause a mild respiratory illness and are not thought to cause epidemics. Influenza A viruses are further categorized into subtypes on t
Keywords	Flu; Matrix protein M1;