



Magic™ Anti-Reovirus Polyclonal antibody (DPATB-H83532)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Reovirus Polyclonal Antibody
Target	Reovirus
Immunogen	Reovirus
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Reovirus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	500 µl
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Reovirus is an acronym for Respiratory Enteric Orphan Virus. The reovirus is a double stranded

RNA virus and is non-pathogenic naturally occurring virus, and not usually associated with any particular illness. Almost 100 per cent of the population show si

Keywords

Respiratory Enteric Orphan Virus,
