



Anti-RSV Polyclonal antibody (DPATB-H83535)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Respiratory Syncytial Virus Polyclonal Antibody
Target	RSV
Immunogen	Viral lysate of Human RSV isolate.
Isotype	IgG
Source/Host	Goat
Species Reactivity	RSV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB, ELISA, Conjugation
Cellular Localization	Virion. Host cytoplasm
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en

Keywords

RS virus,RSV,
