



Anti-SARS-CoV Polyclonal antibody (DPATB-H83537)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-SARS leader Polyclonal Antibody
Target	SARS-CoV
Immunogen	Synthetic peptide: TNHGHKVVVEL VAEMDGIQYG RSGITLGVLV PHVGETPIAY RNVLLRKNGN KGAGGHSYGI DLKSYDLGDE LGTDPIEDYE QNWNTKHGSG ALRELTRELN , corresponding to amino acids 79-180 of SARS TNHGHKVVVELVAEMDGIQYGRSGITLGVLVPHVGETPIAY RNVLLRKNGNKGAGGHSYGIDLKSYDLGDELGTDPIDY
Isotype	IgG
Source/Host	Mouse
Species Reactivity	SARS-CoV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	SARS is not claimed to have been eradicated (unlike smallpox), as it may still be present in its natural host reservoirs (animal populations) and may return to the human population.
Keywords	Human SARS coronavirus,ORF1ab polyprotein,Replicase polyprotein 1ab,