



Magic™ Anti-Sendai virus Polyclonal antibody (DPATB-H83546)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Sendai virus Polyclonal Antibody
Target	Sendai virus
Immunogen	Sendai virus
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Sendai virus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-Fr, ELISA, Neut
Format	Liquid
Size	1 ml
Buffer	40% Glycerol, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Sendai virus is an RNA paramyxovirus parainfluenza type 1 virus that is primarily a disease of young mice; however, it can and does affect several species of small rodents. Sendai virus is immunosuppressive and can have an immediate, as well as a long ter

Keywords

SV,
