



Anti-WNV Core Protein Polyclonal antibody (DPATB-H83556)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-West Nile Virus Core Polyclonal Antibody
Target	WNV Core Protein
Immunogen	Fusion protein: MSKKPGGPGKSRVNMMLKRG, corresponding to amino acids 1/20 of West Nile Virus. MSKKPGGPGKSRVNMMLKRG Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Mouse
Species Reactivity	WNV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The virion of this ssRNA + strand virus is a nucleocapsid covered by a lipoprotein envelope. The nucleocapsid is a complex of Capsid protein C and mRNA.
Keywords	Capsid protein C,