



Anti-WNV PreM Polyclonal antibody (DPATB-H83563)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-West Nile Virus preM Polyclonal Antibody
Target	WNV PreM
Immunogen	Synthetic peptide corresponding to 15 amino acids near the middle of the West Nile Virus Matrix precursor protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	WNV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, ICC/IF, FuncS
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

West Nile Virus is part of the flavivirus family, many of which are capable of causing serious illness in humans. The ssRNA + strand virus is enveloped and spherical with a diameter of 40-50 nm, similar to all flaviviruses. The virion of this virus is a n

Keywords

pre M protein, West Nile virus, WNVgp1,
