



Anti-HSV Polyclonal antibody (DPATB-H83569)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Herpes simplex virus Polyclonal Antibody
Target	HSV
Immunogen	Herpes simplex virus Strain F.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC, IF
Format	Liquid
Size	1 ml
Buffer	pH: 7,20Preservative: 0.1% Sodium azideConstituent: 99% PBS
Preservative	0.1% Sodium Azide
Storage	Store at 4°C or at -20°C for long term storage.

BACKGROUND

Introduction	There are two main types of Herpes Simplex Virus (HSV), 1 and 2. The primary difference between the two types is in where they typically establish latency in the body. HSV1 usually
---------------------	--

establishes latency in the trigeminal ganglion, a collection of nerve cell

Keywords

HSV;
