



Anti-CrPV Polyclonal antibody (DPATB-H83574)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Cricket paralysis virus Polyclonal Antibody
Target	CrPV
Immunogen	Synthetic peptide corresponding to C terminal residues of Cricket paralysis virus
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	CrPV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	50% Glycerol, PBS
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

BACKGROUND

Introduction

Cricket paralysis virus (CrPVgp1) belongs to the picornaviruses polyprotein family.

Keywords

CrPVgp2; CrPVgp2 structural polyprotein; Genome polyprotein; Structural polyprotein;
