



Anti-Simian virus VP1 Polyclonal antibody (DPATB-H83576)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-SV40 VP1 Polyclonal Antibody
Target	Simian virus VP1
Immunogen	Recombinant fusion protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Simian virus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IP, ICC/IF
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

VP1 is one of three simian virus 40 (SV40) capsid proteins. VP1 assembles into pentamers which have a central cavity that contains a copy of one of the other capsid proteins, VP2 or VP3. This complex is then imported into the cell nucleus. Here, 72 of the

Keywords

Capsid protein VP1,
