



Anti-HCV Nonstructural Protein 3, Nonstructural Protein 4 Polyclonal antibody (DPATB-H83600)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Hepatitis C Virus Core NS3 + NS4 Polyclonal Antibody
Target	HCV Nonstructural Protein 3, Nonstructural Protein 4
Immunogen	Recombinant fragment (polypeptide containing core and NS3 and NS4 regions).
Isotype	IgG
Source/Host	Goat
Species Reactivity	HCV
Conjugate	Unconjugated
Applications	Conjugation, ELISA
Cellular Localization	Endoplasmic reticulum and Mitochondrial
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	HCV is a positive, single stranded RNA virus in the Flaviviridae family. The genome is approximately 10,000 nucleotides and encodes a single polyprotein of about 3,000 amino acids. The polyprotein is processed by host cell and viral proteases into three m
Keywords	HCV core and non-structural regions 3 + 4; helicase NS3; Hepatitis C Virus; Non structural protein 4A; Non structural protein 4B; NS3P; NS4A; NS4B; NTPase; p27; p70; p8; Serine protease;