



# Anti-PIV type 3 Polyclonal antibody (DPATB-H83608)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Parainfluenza Virus type 3 Polyclonal Antibody
<b>Target</b>	PIV type 3
<b>Immunogen</b>	Human isolate, type3
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	PIV
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC, IF
<b>Positive Control</b>	Parainfluenza 3 infected cells.
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Infectious diseases are the leading cause of death worldwide. AIDS, tuberculosis (TB), malaria, diarrhoeal and respiratory infections account for 82% of deaths caused by infectious disease. As many infectious diseases are controlled, new diseases emerge a

---

**Keywords**

HPIV3; Parainfluenza 3 virus; PIV3;

---