



Magic™ Anti-Rotavirus Polyclonal antibody (DPATB-H83613)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Rotavirus Polyclonal Antibody
Target	Rotavirus
Immunogen	Rotavirus.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Rotavirus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	500 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Rotaviruses, members of the family Reoviridae, are a major cause of diarrhoea in young mammals. Rotavirus infections also result in economic losses in agriculture due to diarrhoea in calf, pig, sheep, and poultry rearing. Diarrhoea (or scours) due to the

Keywords

Major inner capsid protein VP6,VP6,
