



Magic™ Chicken Anti-Protein L polyclonal antibody [FITC] (DPATB-H83614)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Protein L Polyclonal Antibody, FITC
Target	Protein L
Immunogen	Purified recombinant bacterial Protein L.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Bacteria
Purification	Protein L purified
Conjugate	FITC
Applications	WB, ELISA
Cellular Localization	Cell surface
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.2
Preservative	0.05% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Protein L is a 36,000 dalton immunoglobulin-binding protein isolated from the bacteria *Peptostreptococcus magnus*. Unlike Protein A and Protein G, which bind to the Fc region of immunoglobulins (antibodies), Protein L binds antibodies through light chain i

Keywords

ProteinL; RPL;
