



Anti-*P. aureofaciens* Polyclonal antibody (DPATB-H83615)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat <i>Pseudomonas aureofaciens</i> (conjugated) Polyclonal Antibody
Target	<i>P. aureofaciens</i>
Immunogen	<i>Pseudomonas aureofaciens</i> total proteins
Isotype	IgG
Source/Host	Rat
Species Reactivity	<i>P. aureofaciens</i>
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA, WB
Positive Control	<i>Pseudomonas aureofaciens</i> bacteria
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH ₂ O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Ps. aureofaciens is a naturally occurring bacterium that is widespread in soil and found near plant roots. The bacterium controls certain fungi that cause disease in plants. It is used for controlling various fungal diseases of golf course turf. Large amo

Keywords

Pseudomonas aureofaciens;
