



Anti-P. aureofasciens Polyclonal antibody (DPATB-H83616)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Pseudomonas aureofasciens Polyclonal Antibody
Target	P. aureofasciens
Immunogen	Pseudomonas aureofaciens total proteins
Isotype	IgG
Source/Host	Rat
Species Reactivity	P. aureofasciens
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA, WB
Positive Control	Psuedomonas aureofaciens bacteria
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Ps. aureofaciens is a naturally occurring bacterium that is widespread in soil and found near plant roots. The bacterium controls certain fungi that cause disease in plants. It is used for controlling various fungal diseases of golf course turf. Large amo
Keywords	Pseudomonas aureofasciens,