



# Magic™ Chicken Anti-Protein L Polyclonal antibody [Biotin] (DPATB-H8361B)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant Protein L
<b>Isotype</b>	IgY
<b>Source/Host</b>	Chicken
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Affinity purified.
<b>Conjugate</b>	Biotin
<b>Applications</b>	ELISA, WB Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Reconstitution</b>	Rehydrate by adding 1.0 mL of reagent quality water to the product vial and add preservative.
<b>Format</b>	Lyophilized
<b>Size</b>	1 mg
<b>Buffer</b>	PBS
<b>Preservative</b>	None
<b>Storage</b>	Lyophilized material should be kept at 4°C until rehydrated and will remain stable for a minimum of 2 years. Liquid solution stored at 2-8°C is stable for a minimum of 1 year.

Ship

Wet ice

---