



Anti-17 alpha-Hydroxyprogesterone polyclonal antibody (DPBT-65716SH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti Human 17-Alpha-HydroxyProgesterone-3-CMO Sheep Anti Human 17-Alpha-HydroxyProgesterone-3-CMO
Immunogen	17-Alpha-Hydroxyprogesterone-3-CMO conjugated to KLH
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC, ELISA, WB
Format	Serum - liquid
Size	1 ml
Preservative	0.09% Sodium Azide
Storage	Store at +4 centigrade or at -20 centigrade if preferred. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

BACKGROUND

Introduction	17-Alpha-Hydroxyprogesterone-3-CMO, a progestogen which may reduce the risk of preterm
---------------------	--

delivery in pregnant women.

Keywords

17-ALPHA-HYDROXYPROGESTERONE-3-CMO;
