



Anti-Alpha Neoendorphin polyclonal antibody (DPBT-66731RA)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti Alpha NEO-ENDORPHIN Rabbit Anti Alpha NEO-ENDORPHIN
Immunogen	Alpha neo-Endorphin: Tyr-Gly-Gly-Phe-Leu-Arg-Lys-Tyr-Pro-Lys (porcine).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Pig
Conjugate	Unconjugated
Applications	IHC, ELISA, FC, IP, RIA, WB
Format	Serum - liquid
Size	100 µl
Buffer	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4 °C or at -20 °C if preferred. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

BACKGROUND

Introduction	Neoendorphins are a class of naturally-occurring opioid peptides derived from the proteolytic cleavage of prodynorphin. They include α -neoendorphin and β -neoendorphin.
---------------------	--

Keywords

Neoendorphin; Neo-endorphin
