



# Rabbit Anti-HBsAg (subtypes ad & ay) Polyclonal antibody (DPBT-66855RH)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti Hepatitis B Surface AntigenRabbit Anti Hepatitis B Surface Antigen
<b>Immunogen</b>	Hepatitis B subtypes ad and ay.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Virus
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, FC, IP, IHC, WB
<b>Format</b>	Purified IgG - liquid
<b>Concentration</b>	IgG concentration 4.0 mg/ml
<b>Size</b>	1 ml
<b>Buffer</b>	Phosphate buffered saline
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4 °C or at -20 °C if preferred. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

## BACKGROUND

**Introduction**

Hepatitis B virus, abbreviated HBV, is a species of the genus Orthohepadnavirus, which is likewise a part of the Hepadnaviridae family of viruses. This virus causes the disease hepatitis B.

**Keywords**

HBV pre S2 Antigen; HBV pre S2 protein; Hepatitis B Virus pre S2 protein; PreS2/S antigen; S; Hepatitis B Presurface Antigen; Hepatitis B Presurface; Hepadnaviridae; Orthohepadnavirus; Hepatitis B virus; HBV