



Metformin HCl Standard solution (DWT108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Metformin Standard
-------------------------	---

Conjugate	N/A
------------------	-----

Cas.No	1115-70-4
---------------	-----------

Molecular Weight	165.62g/mol
-------------------------	-------------

Format	Liquid
---------------	--------

Concentration	1.0 mg/mL
----------------------	-----------

Size	1 ml
-------------	------

Buffer	Methanol
---------------	----------

Storage	Store unopened in freezer(-10 °C to 25 °C)
----------------	--

BACKGROUND

Introduction	Metformin is an agent belonging to the biguanide class of antidiabetics with antihyperglycemic activity. It is mainly used for the treatment of type 2 diabetes and can be used alone or in combination with sulfonylureas, thiazolidinediones or other hypoglycemic agents. Up to now, metformin is the most widely used medication for diabetes taken by mouth.
---------------------	---

Keywords	Metformine; Metformin; Dimethylbiguanid; 1,1-Dimethylbiguanide
-----------------	--