



## Codeine Standard solution (DWT22)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	1.0 mg/mL analytical Codeine Standard
-------------------------	---------------------------------------

<b>Conjugate</b>	N/A
------------------	-----

<b>Cas.No</b>	76-57-3
---------------	---------

<b>Molecular Weight</b>	299.36g/mol
-------------------------	-------------

<b>Format</b>	Liquid
---------------	--------

<b>Concentration</b>	1.0 mg/mL
----------------------	-----------

<b>Size</b>	1 ml
-------------	------

<b>Buffer</b>	Methanol
---------------	----------

<b>Storage</b>	Store unopened in freezer(-10°C to 25°C)
----------------	--

### BACKGROUND

<b>Introduction</b>	Codeine is an opioid analgesic used to treat moderate to severe pain when the use of an opioid is indicated.
---------------------	--

<b>Keywords</b>	Ascomp; Cheratussin; Codeine polistirex; Codeinum; Methylmorphine; Codeine Contin; Covan; Dimetane
-----------------	--