



Clonidine Standard solution (DWT23)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Clonidine Standard
------------------	---

Conjugate	N/A
-----------	-----

Cas.No	4205-90-7
--------	-----------

Molecular Weight	230.09g/mol
------------------	-------------

Format	Liquid
--------	--------

Concentration	1.0 mg/mL
---------------	-----------

Size	1 ml
------	------

Buffer	Methanol
--------	----------

Storage	Store unopened in freezer(-10°C to 25°C)
---------	--

BACKGROUND

Introduction	Clonidine is an alpha 2-adrenoceptor agonist. Clonidine is extensively used as an adjunct to spinal morphine and is suggested to be particularly effective for neuropathic pain; its clinical use is limited by side effects such as sedation and hypotension.
--------------	--

Keywords	Clonidine; alpha 2-adrenoceptor agonist
----------	---