



Meperidine Standard solution (DWT31)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Meperidine Standard
Conjugate	N/A
Cas.No	57-42-1
Molecular Weight	247.33g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Meperidine is an opioid agonist with analgesic and sedative properties used to manage severe pain.
Keywords	Demerol; Meperidine; Pethidine