



## Felbamate Standard solution (DWT43)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Product Overview</b>	2.0 mg/mL analytical Felbamate Standard
-------------------------	---

<b>Conjugate</b>	N/A
------------------	-----

<b>Cas.No</b>	25451-15-4
---------------	------------

<b>Molecular Weight</b>	238.24g/mol
-------------------------	-------------

<b>Format</b>	Liquid
---------------	--------

<b>Concentration</b>	2.0 mg/mL
----------------------	-----------

<b>Size</b>	1 ml
-------------	------

<b>Buffer</b>	Methanol
---------------	----------

<b>Storage</b>	Store unopened in freezer(-10°C to 25°C)
----------------	--

### BACKGROUND

<b>Introduction</b>	Felbamate is an anticonvulsant drug used in the treatment of epilepsy. In particular, in the adult patient population, it can be employed to treat partial seizures (with and without generalization). Alternatively, it is used to treat partial and generalized seizures associated with Lennox-Gastaut syndrome in children. It has a weak inhibitory effect on GABA receptor binding sites.
---------------------	---

<b>Keywords</b>	Felbamate; Felbamato; Felbamatum
-----------------	----------------------------------