



Temazepam Standard solution (DWT53)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Temazepam Standard
Conjugate	N/A
Cas.No	846-50-4
Molecular Weight	300.74g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Temazepam is approved for the short-term treatment of insomnia. In addition, temazepam has anxiolytic (anti-anxiety), anticonvulsant, and skeletal muscle relaxant properties.
Keywords	Temazepam; Restoril; Normison