



Imipramine HCl Standard solution (DWT62)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Imipramine Standard
Conjugate	N/A
Cas.No	113-52-0
Molecular Weight	316.87g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Imipramine (sold as Antidepressin, Deprimin, Deprinol, Depsol, Depsonil, Dynaprin, Eupramin, Imipramil, Irmin, Janimine, Melipramin, Surplix, Tofranil), also known as melipramine, is an antidepressant medication, a tricyclic antidepressant of the dibenzazepine group. Imipramine is mainly used in the treatment of major depression and enuresis (inability to control urination).
Keywords	Antidepressin; Deprimin; Deprinol; Depsol; Depsonil; Dynaprin; Eupramin; Imipramil; Irmin; Janimine; Melipramin; Surplix; Tofranil; melipramine; Imipramine