



Trazodone HCl Standard solution (DWT65)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Trazodone Standard
Conjugate	N/A
Cas.No	25332-39-2
Molecular Weight	408.32g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Trazodone is a serotonin uptake inhibitor used to treat major depressive disorder.
Keywords	Trazodone; Desyrel; Oleptro; Beneficat; Deprax; Desirel