



Amiodarone HCl Standard solution (DWT70)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Amiodarone Standard
Conjugate	N/A
Cas.No	19774-82-4
Molecular Weight	681.77g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Amiodaron is an antianginal and antiarrhythmic drug. It increases the duration of ventricular and atrial muscle action by inhibiting Na,K-activated myocardial adenosine triphosphatase. There is a resulting decrease in heart rate and in vascular resistance.
Keywords	Nexterone; Pacerone; Amiodarone