



(-)-Levamisole HCl Standard solution (DWT77)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	1.0 mg/mL analytical Levamisole Standard
Conjugate	N/A
Cas.No	16595-80-5
Molecular Weight	240.75g/mol
Format	Liquid
Concentration	1.0 mg/mL
Size	1 ml
Buffer	Methanol
Storage	Store unopened in freezer(-10°C to 25°C)

BACKGROUND

Introduction	Levamisole is levamisole is a nicotinic receptor agonist used to treat helminth infections and some skin infections.
Keywords	Levamisole