



(±)-Flecainide Standard solution (DWT96)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview 1.0 mg/mL analytical Flecainide Standard

Conjugate N/A

Cas.No 54143-55-4

Molecular Weight 414.34g/mol

Format Liquid

Concentration 1.0 mg/mL

Size 1 ml

Buffer Methanol

Storage Store unopened in freezer(-10 °C to 25 °C)

BACKGROUND

Introduction Flecainide is a potent anti-arrhythmia agent, effective in a wide range of ventricular and atrial arrhythmias and tachycardias. Paradoxically, however, in myocardial infarct patients with either symptomatic or asymptomatic arrhythmia, flecainide exacerbates the arrhythmia and is not recommended for use in these patients.

Keywords (±)-flecainide; Flecaïne; Flecainide; N-(2-Piperidinylmethyl)-2,5-bis(2,2,2-trifluoroethoxy)benzamide